

Your expected body is possible,
get it now !

CryoSlim[®] Hybrid

Cryolipolysis and shock waves




CryoSlim Hybrid® : the Cryolipolysis with guarantee of results

The CryoSlim Hybrid® is the unique slimming device that guarantees the results for patients.

Intuitive and interactive, the CryoSlim Hybrid® communicates the expected results to the practitioner within the 10 first seconds of treatment.


Managed by the Artificial Intelligence USER ASSIST®, the CryoSlim Hybrid® analyses the fat treated and optimizes automatically the session parameters according to the morphology of each area.




Cryolipolysis : Fat elimination by cold

Cryolipolysis is currently the best technique for fat reduction with durable effects. Alternative to the liposuction, the Cryolipolysis works following a natural process which targets the localized adipocytes cells (saddle bags, abdomen, back, legs..).


Adipocytes are sensitive to the change of temperature and especially to cold: an intensive cold exposure generates the progressive elimination of those fat cells.



Before



Treatment



After

Cryolipolysis new generation, the CryoSlim Hybrid® treats the fat cells by combined cold, hot and vacuum actions.

USER ASSIST® : Artificial Intelligence


Thanks to its interactive abilities USER ASSIST® guarantees the quality of the expected results.

USER ASSIST® is an Artificial Intelligence which reacts in live with the patient's body. Fat volume, fat density and fat composition are analyzed by USER ASSIST®.

A cross referring data adjusts automatically the Cryolipolysis parameters to be perfectly suitable for each treatment areas and patient's morphologies.



USER ASSIST® : Diagnose and communicate with the practitioner




USER ASSIST® is the only one assistance technology ables to guarantee the results for patients and actively works for their safety.

Ovoid Cryodes® : Combine efficacy and safety



Morphological with its ovoid shape, the Cryodes® of the CryoSlim Hybrid® are exclusive and unique in the world. Striated over its contact surface they ensure an homogeneous distribution of suction for optimal comfort.



Pictures made by thermal camera FLIR-360 T high resolution immediately after one Cryolipolysis session done with the CryoSlim Hybrid®.

Patented for their cooling system, the Cryodes® manages two times more thermal extractors than all other Cryolipolysis devices.

This technical exclusivity concentrates the efficacy at treatment temperatures from 0°C to -5°C. Healthy for the organism those temperatures preserve the body from uncomfot post Cryolipolysis.

This is evidenced by the 93% of patients satisfaction*.

CryoSlim Hybrid® technologies



Shock Waves®

The Shock Waves® is used in order to break the fat fibrosis and make the fat more flexible.

Fat enter easier inside the Cryode® and generates better results for all patients and especially patients with low level of Cryolipolysis results.

Indications : before and/or after the Cryolipolysis



Large Cryode®

Indications : abdomen, saddle bags, hips, back



Small Cryode®

Indications : arms, knees, legs

Clinical studies done with Cryolipolysis devices demonstrate an average fat elimination from 15 % to 25 % of fat per session.

Temperatures used are mainly between -7°C to -12°C .

Synthesis of 10 years of experience in Cryolipolysis devices, the CryoSlim Hybrid reaches better results from 0°C treatment temperature :


- From 29% to 39% of fat volume removed per session*
- 93% of patients are satisfied*
- 90% of patients recommend the session*
- Painless*

First results appear 3 weeks after the Cryolipolysis session.


The slimming effects continue during 2 to 3 months post session.

* Adjadj L, SidAhmed-Mezi M, Mondoloni M, Meningaud JP, Hersant B. Assessment of the Efficacy of Cryolipolysis on Saddlebags: A Prospective Study of 53 Patients. *Plast Reconstr Surg.* 2017 Jul;140(1):50-57

al results



Before



After



Before



After

Courtesy Dr. Ariana

Results after 1 session of CryoSlim Hybrid®




Before



After



Before



After

Courtesy Unique Skin

Results after 3 session of CryoSlim Hybrid®

Non invasive and without any pain, patients maintain their usual way of life, professional and social activities.



+ 33 4 66 32 31 16 • aamsworld.fr • info@aamsworld.com

Cryode®	Ovoid / 4 Peltiers (heat extractor)
Vacuum	Automatically managed by USER ASSIST from 0 to 600 mbar
Temperature	+ 37 °C / - 5° C
Screen	High resolution 7 inch tactile screen
Tension	Electrical tension 280 W / Liquid cooling system
Sizes	60cm x 35cm x 118cm

Designed and made in our factories in France, the CryoSlim Hybrid® is followed and certified by the GMED

